

Product Datasheet

VWF antibody (orb536997)

Description	Rabbit polyclonal antibody to VWF
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human von Willebrand factor (VWF)(NP_000543.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	309KDa
Uniprot ID	P04275
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.