

Product Datasheet

CTNNB1 (Phospho-S675) antibody (orb537006)

Description	Rabbit polyclonal antibody to CTNNB1 (Phospho-S675)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S675 of human beta Catenin(NP_001091679.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	100kDa
Uniprot ID	P35222
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.

