


Product Datasheet

GRIN1 (Phospho-S890) antibody (orb537008)

Description	Rabbit polyclonal antibody to GRIN1 (Phospho-S890)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S890 of human NMDAR1(NP_015566.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	140kDa
Uniprot ID	Q05586
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of C6 ...