



Product Datasheet

Histone HIST1H3A (Phospho-T3) antibody (orb537009)

Description	Rabbit polyclonal antibody to Histone HIST1H3A (Phospho-T3)
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T3 of human HIST1H3A(NP_003520.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	17KDa
Uniprot ID	P68431
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Western blot analysis of extracts of C6 ...