


Product Datasheet

PERK (Phospho-T982) antibody (orb537011)

Description	Rabbit polyclonal antibody to PERK (Phospho-T982)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding T982 of human PERK.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	170KDa
Uniprot ID	Q9NZJ5
Entrez	9451
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of C6 ...