


Product Datasheet

BACE2 antibody (orb537014)

Description	Rabbit polyclonal antibody to BACE2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human BACE2(NP_036237.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56kDa
Uniprot ID	Q9Y5Z0
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of DU145 cells.