


Product Datasheet

CDK4 (Phospho-T172) antibody (orb537032)

Description	Rabbit polyclonal antibody to CDK4 (Phospho-T172)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T172 of human CDK4(NP_000066.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	33KDa
Uniprot ID	P11802
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HCT...