



Product Datasheet

CYP19A1 antibody (orb537052)


Description	Rabbit polyclonal antibody to CYP19A1
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Aromatase (CYP19A1)(NP_000094.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48kDa
Uniprot ID	P11511
Dilution Range	WB: 1:1000-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat ovar...



Immunohistochemical staining of human pl...



Immunohistochemical staining of human co...