


Product Datasheet

FHOD1 antibody (orb537058)

Description	Rabbit polyclonal antibody to FHOD1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human FHOD1(NP_037373.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	145kDa
Uniprot ID	Q9Y613
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa...