


Product Datasheet

IAPP antibody (orb537068)

Description	Rabbit polyclonal antibody to IAPP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	A synthetic peptide of human IAPP(NP_000406.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	9kDa
Uniprot ID	P10997
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa...