


Product Datasheet

Loxl4 antibody (orb53707)

Description	Loxl4 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Lysyl oxidase homolog 4 protein (26-757AA)
Target	Loxl4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q924C6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat brain tissue...