



Product Datasheet

IGF2R antibody (orb537071)


Description	Rabbit polyclonal antibody to IGF2R
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Cation-independent M6PR (Cation-independent M6PR (IGF2R))(NP_000867.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	274KDa
Uniprot ID	P11717
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of C6 cells ...



Immunohistochemical staining of human ga...