



Product Datasheet

KRT7 antibody (orb537077)


Description	Rabbit polyclonal antibody to KRT7
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Cytokeratin 7 (Cytokeratin 7 (KRT7))(NP_005547.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51kDa
Uniprot ID	P08729
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse bronchus



Immunofluorescence analysis of rat bronchus



Immunofluorescence analysis of human lung