


Product Datasheet

MAP2K1 (Phospho-pS221) antibody (orb537108)

Description	Rabbit polyclonal antibody to MAP2K1 (Phospho-pS221)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S221 of human MEK1(NP_002746.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	43KDa
Uniprot ID	Q02750
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa cells.