



Product Datasheet

HSD11B2 antibody (orb537180)


Description	Rabbit polyclonal antibody to HSD11B2
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human HSD11B2(NP_000187.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Uniprot ID	P80365
Dilution Range	WB: 1:200-1:1000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of human st...



Immunohistochemical staining of human li...