


Product Datasheet

CLINT1 antibody (orb537206)

Description	Rabbit polyclonal antibody to CLINT1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic Peptide of human CLINT1
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	64kDa
Uniprot ID	Q14677
Entrez	9685
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HL-...