



Product Datasheet

NTF4 antibody (orb537217)


Description	Rabbit polyclonal antibody to NTF4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NTF4 (NP_006170.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	35kDa
Uniprot ID	P34130
Entrez	4909
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human prostate tissue.



Immunohistochemical staining of human colon tissue.



Western blot analysis of extracts of HT-1080 cells.