

Product Datasheet

EME2 antibody (orb537226)

Description	Rabbit polyclonal antibody to EME2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human EME2(NP_001244299.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41kDa
Uniprot ID	A4GXA9
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.