


Product Datasheet

Jak1 (Phospho-Y1022/1023) antibody (orb537234)

Description	Rabbit polyclonal antibody to Jak1 (Phospho-Y1022/1023)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic phosphorylated peptide around Y1022 & Y1023 of human JAK1(NP_002218.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	130KDa
Uniprot ID	P23458
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HT-29