



Product Datasheet

TP53 (Phospho-S33) antibody (orb537236)


Description	Rabbit polyclonal antibody to TP53 (Phospho-S33)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthetic phosphorylated peptide around S33 of human p53(NP_000537.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	55kDa
Uniprot ID	P04637
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat ovar...



Immunohistochemical staining of rat test...



Immunohistochemical staining of rat live...