



Product Datasheet

TNFRSF6B antibody (orb537245)


Description	Rabbit polyclonal antibody to TNFRSF6B
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TNFRSF6B(NP_003814.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33KDa
Uniprot ID	O95407
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat kidney tissue.



Immunohistochemical staining of human kidney tissue.



Immunohistochemical staining of human tissue.