


Product Datasheet

MMP12 antibody (orb537256)

Description	Rabbit polyclonal antibody to MMP12
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human MMP12(NP_002417.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45KDa
Uniprot ID	P39900
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of hum...