



Product Datasheet

AZGP1 antibody (orb537260)


Description	Rabbit polyclonal antibody to AZGP1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Zinc- α 2-glycoprotein (ZAG/Zinc- α 2-glycoprotein (ZAG/AZGP1))(NP_001176.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41kDa
Uniprot ID	P25311
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lung tissue.



Immunohistochemical staining of human prostate tissue.



Immunohistochemical staining of human colon tissue.