



Product Datasheet

PPBP antibody (orb537262)

Description	Rabbit polyclonal antibody to PPBP
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PPBP(NP_002695.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	P02775
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human li...



Western blot analysis of extracts of hum...