


Product Datasheet

Jak2 (Phospho-Y1008) antibody (orb537295)

Description	Rabbit polyclonal antibody to Jak2 (Phospho-Y1008)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y1008 of human Phospho-JAK2-Y1008(NP_004963.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	120kDa
Uniprot ID	O60674
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Jur...