



Product Datasheet

MAPK3 (Phospho-T202) antibody (orb537297)


Description	Rabbit polyclonal antibody to MAPK3 (Phospho-T202)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic phosphorylated peptide around T202 of human ERK1/ERK2(NP_002737.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	42KDa/44KDa
Uniprot ID	P27361 , P28482
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:100-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human br...



Western blot analysis of extracts of Jur...