

Product Datasheet

MAPK3 (Phospho-Y204) antibody (orb537298)

Description Rabbit polyclonal antibody to MAPK3 (Phospho-Y204)

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, IP, WB

Immunogen A synthetic phosphorylated peptide around Y204 of human ERK1(NP_002737.2).

Preservatives PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only


Isotype IgG

MW 42KDa/44KDa


Uniprot ID **P27361**

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:100-1:200, IP: 1:50-1:100


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of rat ovar...



Immunohistochemical staining of rat test...



Immunohistochemical staining of rat live...