



Product Datasheet

PLCG1 (Phospho-S1248) antibody (orb537300)

Description	Rabbit polyclonal antibody to PLCG1 (Phospho-S1248)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around S1248 of human PLC gamma 1 (PLCG1) (NP_002651.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	149kDa
Uniprot ID	P19174
Entrez	5335
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:10-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...



Western blot analysis of extracts of Jur...