



Product Datasheet

STMN1 (Phospho-S25) antibody (orb537304)

Description	Rabbit polyclonal antibody to STMN1 (Phospho-S25)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around S25 of human Stathmin 1(NP_005554.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	17kDa
Uniprot ID	P16949
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Western blot analysis of extracts of Jurkat cells.