

Product Datasheet

TBC1D1 antibody (orb537311)

Description	Rabbit polyclonal antibody to TBC1D1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human TBC1D1(NP_055988.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150KDa
Uniprot ID	Q86T10
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.