


Product Datasheet

STAT5A (Phospho-Y694) antibody (orb537332)

Description	Rabbit polyclonal antibody to STAT5A (Phospho-Y694)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y694 of human STAT5A(NP_001275647.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	90kDa
Uniprot ID	P42229
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of K56...