

Product Datasheet

STMN1 (Phospho-S38) antibody (orb537334)

Description Rabbit polyclonal antibody to STMN1 (Phospho-S38)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen A phospho specific peptide corresponding to residues surrounding S38 of human Stathmin 1

Preservatives Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only

Isotype IgG

Clonality Polyclonal

Purity Affinity purification


MW 19KDa

Uniprot ID **P16949**

Entrez **3925**

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of K-5...