



Product Datasheet

ATP1A antibody (orb537340)


Description	Rabbit polyclonal antibody to ATP1A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ATP1A (NP_000692.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	113kDa
Uniprot ID	P05023
Entrez	476477478480
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human st...



Immunohistochemical staining of human pr...



Western blot analysis of extracts of LO2...