


Product Datasheet

PRG4 antibody (orb537358)

Description	Rabbit polyclonal antibody to PRG4
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PRG4(NP_001121180.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150kDa
Uniprot ID	Q92954
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of LO2...