



Product Datasheet

POLR2A (Phospho-S5) antibody (orb537385)


Description	Rabbit polyclonal antibody to POLR2A (Phospho-S5)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S5 of human POLR2A CTD repeatYSPTSPS.
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	270kDa
Uniprot ID	P24928
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF7 cell...



Immunohistochemical staining of mouse li...



Immunohistochemical staining of human CO...