


Product Datasheet

CAPNS2 antibody (orb537415)

Description	Rabbit polyclonal antibody to CAPNS2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CAPNS2 (NP_115706.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	32kDa
Uniprot ID	Q96L46
Entrez	84290
Dilution Range	WB: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...