



Product Datasheet

GRIN2B antibody (orb537437)


Description	Rabbit polyclonal antibody to GRIN2B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human GRIN2B(NP_000825.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	166KDa
Uniprot ID	Q13224
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse brain...



Immunofluorescence analysis of rat brain...



Western blot analysis of extracts of mou...