

Product Datasheet

PIK3C3 antibody (orb53747)

Description	PIK3C3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Phosphatidylinositol 3-kinase catalytic subunit type 3 protein (700-850AA)
Target	PIK3C3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NEB9
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.