


Product Datasheet

ARMS2 antibody (orb537506)

Description	Rabbit polyclonal antibody to ARMS2
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ARMS2(NP_001093137.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11kDa
Uniprot ID	P0C7Q2
Dilution Range	WB: 1:200-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mouse eye.