


Product Datasheet

TFAP2A antibody (orb537518)

Description	Rabbit polyclonal antibody to TFAP2A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Activator protein 2 (AP-2/Activator protein 2 (AP-2/TFAP2A)) (NP_001027451.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	50KDa
Uniprot ID	P05549
Entrez	7020
Dilution Range	WB: 1:200-1:1000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...