


Product Datasheet

CAMK2D antibody (orb537527)

Description	Rabbit polyclonal antibody to CAMK2D
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human CaMKII delta (NP_742125.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	50kDa
Uniprot ID	Q13557
Entrez	817
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mouse heart tissue.