



Product Datasheet

METTL20 antibody (orb537546)

Description	Rabbit polyclonal antibody to METTL20
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human METTL20(NP_001129335.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33kDa
Uniprot ID	Q8IXQ9
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...



Western blot analysis of extracts of mouse...