


Product Datasheet

HEMK1 antibody (orb537587)

Description	Rabbit polyclonal antibody to HEMK1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human HEMK1(NP_057257.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38kDa
Uniprot ID	Q9Y5R4
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...