



Product Datasheet

IFNL3 antibody (orb537663)


Description	Rabbit polyclonal antibody to IFNL3
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human IFNL3(NP_742151.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24KDa
Uniprot ID	Q8IZI9
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse liver tissue.



Immunohistochemical staining of human tissue.



Western blot analysis of extracts of mouse liver.