


Product Datasheet

KRT4 antibody (orb537666)

Description	Rabbit polyclonal antibody to KRT4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human Cytokeratin 4 (Cytokeratin 4 (KRT4))
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	55kDa
Uniprot ID	P19013
Entrez	3851
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...