



Product Datasheet

LAMP1 antibody (orb537668)


Description	Rabbit polyclonal antibody to LAMP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human LAMP1(NP_005552.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa/120kDa
Uniprot ID	P11279
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 cells...



Immunofluorescence analysis of HeLa cells...



Western blot analysis of extracts of mouse liver...