


Product Datasheet

PER2 antibody (orb537681)

Description	Rabbit polyclonal antibody to PER2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IHC, WB
Immunogen	Recombinant fusion protein of human PER2(NP_073728.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	136KDa
Uniprot ID	O15055
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mouse liver.