


Product Datasheet

CTNNB1 antibody (orb537715)

Description	Rabbit polyclonal antibody to CTNNB1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human beta Catenin(NP_001895.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100KDa
Uniprot ID	P35222
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:50-1:100, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...