


Product Datasheet

LYPLAL1 antibody (orb537761)

Description	Rabbit polyclonal antibody to LYPLAL1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human LYPLAL1(NP_620149.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28kDa
Uniprot ID	Q5VWZ2
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Mou...