


Product Datasheet

UCP2 antibody (orb537817)

Description	Rabbit polyclonal antibody to UCP2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human UCP2(NP_003346.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32KDa
Uniprot ID	P55851
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...